

NEW!

Terumo SurSaver® 23G Mixing Syringe Tray

Maximize the Mix Minimize the Effort

- Easy to draw up viscous allergenic extracts
- Only 0.0021 mL dead space¹
- 98% less antigen waste (see back)
- Smooth plunger
- Clear barrel
- Easy to read scale markings
- Ultra-sharp needle

¹ Testing demonstrated an average of 0.0021 dead space.



Mixing Syringes/Safety and Conventional Allergy Syringes

Ordering Information:

NEW

SurSaver® Mixing Syringe Trays

Packed: 25/tray, 40 trays/case (1,000)

Product Code	Description
SS01A2313T	1cc with 23G x 1/2" regular bevel needle, tray

SurSaver® Mixing Syringes Individually Packed

Packed: 100/box, 5 boxes/case (500)

Product Code	Description
SS01D2313	1cc with 23G x 1/2" regular bevel needle, individually packed



LATEX FREE



LATEX FREE

Comparative Cost Savings Using SurSaver® Mixing Syringes

Antigen Cost per mL	Average Dead Space per Use		Savings per Use*
	Terumo® SurSaver® Mixing Syringe	Terumo® 1cc syringe with a luer slip 23G needle	
	0.0021 mL	0.1030 mL	
\$3.00	\$0.006	\$0.309	\$0.30
\$6.00	\$0.013	\$0.618	\$0.61
\$9.00	\$0.019	\$0.927	\$0.91
\$12.00	\$0.025	\$1.236	\$1.21

* Savings based upon average usage, individual results may vary

Terumo SurSaver® Allergy Mixing Syringes will minimize your antigen waste with an average of 98% less dead space than a similar Terumo slip tip product.

Need Samples?

We've now made it faster and easier for you to obtain samples and literature of our products through your local TMP rep. Just log-on to www.terumoTMP.com and select "Find Your Rep" from our homepage.

Ⓟ and TERUMO are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. SURSAVER is a registered trademark of Terumo Medical Corporation, Somerset, NJ.
© 2008 Rev. 9/08

TERUMO®

TERUMO MEDICAL PRODUCTS